



ISO 9001:2015 Certified Institution

Balaji Institute of Technology & Science

Estd.:2001

Laknepally (V), Narsampet (M), Warangal District - 506 331, Telangana State, India

(AUTONOMOUS)

Accredited by NBA (UG - CE, ME, ECE & CSE) **& NAAC A+ Grade**

(Affiliated to JNT University, Hyderabad and Approved by AICTE, New Delhi)

www.bitswgl.ac.in, email: principal@bitswgl.ac.in, Ph:98660 50044, Fax: 08718-230521

ENGINEERINGCHEMISTRY LAB MANUAL

Course Code : **22CH102BS**
Regulations : **R22**
Class : **B.Tech**
Branch : **CE/EEE/ME/ECE/CSE/CSW/CSM/CSO**

Preparedby

Rizwana Farheen

Assistant Professor in Chemistry,
DepartmentofHumanitiesandSciences



ISO 9001:2015 Certified Institution

Balaji Institute of Technology & Science

Estd.:2001

Laknepally (V), Narsampet (M), Warangal District - 506 331, Telangana State, India

(AUTONOMOUS)

Accredited by NBA (UG - CE, ME, ECE & CSE) & NAAC A+ Grade

(Affiliated to JNT University, Hyderabad and Approved by AICTE, New Delhi)

www.bitswgl.ac.in, email: principal@bitswgl.ac.in, Ph:98660 50044, Fax: 08718-230521

Vision of Department

To establish a centre of excellence in basic sciences such as Mathematics, Physics, Chemistry & Environmental Sciences that provides foundation for engineering studies and also in English Language Communication Skills that helps students to express themselves effectively and to create engineers with proficiency in engineering fundamentals – experimental, analytical, computational and designing abilities.

Mission of Department

- M1: To create academic excellence in fundamental sciences and communication skills for the Students.
- M2: To encourage advanced teaching learning process, quality based knowledge and Quality research at individual, department and institutional level.
- M3: To impart personality development skills to students that will help them to succeed and lead.

Course Objectives

- To bring adaptability to new developments in Engineering Chemistry and to acquire the skills required to become a perfect engineer.
- To include the importance of water in industrial usage, fundamental aspects of battery chemistry, significance of corrosion its control to protect the structures.
- To imbibe the basic concepts of petroleum and its products
- To acquire required knowledge about engineering materials like cement, smart materials and Lubricants.



ISO 9001:2015 Certified Institution

Balaji Institute of Technology & Science

Estd.:2001

Laknepally (V), Narsampet (M), Warangal District - 506 331, Telangana State, India

(AUTONOMOUS)

Accredited by NBA (UG - CE, ME, ECE & CSE) & NAAC A+ Grade

(Affiliated to JNT University, Hyderabad and Approved by AICTE, New Delhi)

www.bitswgl.ac.in, email: principal@bitswgl.ac.in, Ph:98660 50044, Fax: 08718-230521**Program Outcomes**

PO1:	Engineering knowledge: Apply the knowledge of mathematics, science, engineering fundamentals, and an engineering specialization to the solution of complex engineering problems.
PO2:	Problem analysis: Identify, formulate, review research literature, and analyze complex engineering problems reaching substantiated conclusions using first principles of mathematics, natural sciences, and engineering sciences.
PO3:	Design/development of solutions: Design solutions for complex engineering problems and design system components or processes that meet the specified needs with appropriate consideration for the public health and safety, and the cultural, societal, and environmental considerations.
PO4:	Conduct investigation of complex problems: Use research-based knowledge and research methods including design of experiments, analysis and interpretation of data, and synthesis of the information to provide valid conclusions.
PO5:	Modern tool usage: Create, select, and apply appropriate techniques, resources, and modern engineering and IT tools including prediction and modeling to complex engineering activities with an understanding of the limitations.
PO6:	The engineer and society: Apply reasoning informed by the contextual knowledge to assess societal, health, safety, legal and cultural issues and the consequent responsibilities relevant to the professional engineering practice.
PO7:	Environment and sustainability: Understand the impact of the professional engineering solutions in societal and environmental contexts, and demonstrate the knowledge of, and need for sustainable development.
PO8:	Ethics: Apply ethical principles and commit to professional ethics and responsibilities and norms of the engineering practice.
PO9:	Individual and teamwork: Function effectively as an individual, and as a member or leader in diverse teams, and in multidisciplinary settings.
PO10:	Communication: Communicate effectively on complex engineering activities with the engineering community and with society at large, such as, being able to comprehend and write effective reports and design documentation, make effective presentations, and give and receive clear instructions.
PO11:	Project management and finance: Demonstrate knowledge and understanding of the engineering and management principles and apply these to one's own work, as a member and leader in a team, to manage projects and in multidisciplinary environments.
PO12:	Life-long learning: Recognize the need for, and have the preparation and ability to engage in independent and life-long learning in the broadest context of technological change.



ISO 9001:2015 Certified Institution

Balaji Institute of Technology & Science

Estd.:2001

Laknepally (V), Narsampet (M), Warangal District - 506 331, Telangana State, India

(AUTONOMOUS)

Accredited by NBA (UG - CE, ME, ECE & CSE) & NAAC A+ Grade

(Affiliated to JNT University, Hyderabad and Approved by AICTE, New Delhi)

www.bitswgl.ac.in, email: principal@bitswgl.ac.in, Ph:98660 50044, Fax: 08718-230521

CourseOutcomes

CO1:	Understand Water treatment, specifically hardness of water and purification of water by various methods..
CO2:	Determine the conductometry , potentiometry , and pH metry & the concentrations or equivalence points of acids and bases would generally focus on building the students' understanding of these analytical techniques and their practical applications.
CO3:	Understanding Thermoplastics & the properties and characteristics of thermoplastics , with a focus on Nylon 6 , including its flexibility, mold ability, and recyclability. Students will learn about the characteristics and properties of thermosetting plastics, particularly Bakelite .
CO4:	Estimation of physical properties like viscosity and acid value of oil samples. Students will understand the concept of viscosity as a measure of a fluid's resistance to flow, and how it relates to the internal friction of the oil.
CO5:	Illustrating the construction of fuel cells , solar cells , and smart materials focus on equipping students with both theoretical knowledge and practical skills to understand and design these advanced materials and technologies.

ATTAINMENT OF COURSE OUTCOMES & PROGRAM OUTCOMES

COs	PO1	PO2	PO3	PO4	PO5	PO6	PO7	PO8	PO9	PO10	PO11	PO12
CO1	2	2	-	-	3	--	3	-	3	3	-	-
CO2	2	-	-	-	3	-	-	--	3	3	-	-
CO3	3	-	-	-	3	-	2	-	3	3	-	-
CO4	3	-	-	-	3	-	-	-	3	3	-	-
CO5	3	-	-	-	3	-	2	-	3	3	-	-
Average	2.6	2	-	-	3	-	2.3	-	3	3	-	-

1=Slite(Low)

2=Moderate(Medium)

3=Substantial(High)

**BALAJIGROUP OF INSTITUTION
LAKNEPALLY, NARSAMPET, WARANGAL**



CERTIFICATE

Name of the Laboratory.....

*This is to certify that, this is a bonafied record of laboratory experimental work done by ... bearing Hall
Ticket Number B.Tech/MBA/MCA
..... Semester student in the department
of of this college during the year*

(Signature of the concerned teacher) APPROVED/ NOT

APPROVED

External Examiner

Internal Examiner

Department of

INDEX

S.No	ExperimentName	Page No.	Date of Experiment	Date of Submission	Remarks
1	Determination of total hardness of water by complexometric method using EDTA.				
2	Estimation of the concentration of an acid by Conductometry.				
3	Estimation of the amount of Fe+2 by Potentiometry.				
4	Determination of an acid concentration using pH meter.				
5	Preparation of Bakelite.				
6	Preparation Nylon – 6				
7	Estimation of acid value of given lubricant oil.				
8	Estimation of Viscosity of lubricant oil using Ostwald's Viscometer				
9	Determination of rate of corrosion of mild steel in the presence and absence of inhibitor.				
10	Construction of Fuel cell and its working.				
11	Smart materials for Biomedical applications.				
12	Batteries for electrical vehicles..				
13	Functioning of solar cell and its applications				

GENERAL PRECAUTIONS FOR AVOIDING ACCIDENTS IN A CHEMISTRY LABORATORY

- Use full length protective coat/an apron.
- Wear safety goggles for preventing eye injuries by splashing of chemicals.
- A pair of gloves must be used when handling poisonous or toxic chemicals.
- A rubber bulb should be used for pipetting solutions.
- Before leaving the laboratory, wash and clean every apparatus, wipe the table and keep all the apparatus in its original position.
- Laboratory must be kept clean and tidy at all times.
- Each laboratory must be equipped with first aid box.
- Strong acids used for conducting experiments should be diluted before being poured into the sink.
- When the burner is not being used, lower the flame. As soon as the day's work is over close the gas tap.
- Never handle chemicals with your finger. Always use fresh, clean spatula.
- Never try to lift the glass bottle holding the neck alone.
- Don't eat, drink in the laboratory.
- Do not throw solid wastes or filter paper or other wastes into the sink.

LAB SAFETY RULES

- **ALWAYS** wear your apron or protective clothing when working with chemicals.
- **ALWAYS** tie back loose hair.
- **ALWAYS** wear goggles or safety glasses to prevent getting materials in your eyes.
- **ALWAYS** read the labels and heed all warnings.
- **NEVER** eat, drink or smell the chemicals. Rather carefully "fan" the fumes to your nose.
- **NEVER** look directly into a test tube or flask. Look at the contents from the side.
- **NEVER** play around during experiments.
- **ALWAYS** wash your hands after handling lab materials.

VOLUMETRIC ANALYSIS

BASIC CONCEPTS OF VOLUMETRIC ANALYSIS

Chemical analysis of the compounds is carried out in two ways

1. **Qualitative analysis.**
2. **Quantitative analysis.**

Qualitative analysis shows what element a given contains.

Quantitative analysis determines the quantity of a particular component present in substance. It is carried out in two ways

1. **Gravimetric analysis.**
2. **Volumetric analysis.**

Gravimetric analysis involves the estimation of the amount of a given compound from the results of weighing.

Volumetric analysis is based on measuring the volume of the solution of a substance. Terms involved in volumetric analysis

1. **Titration:** The process of finding out the volume of one of the solutions required to react completely with a definite volume of the other solution of known concentration is called titration.
2. **Titrant:** The solution of known strength is called titrant.
3. **Titrate:** The solution whose concentration is to be estimated.
4. **Indicator:** The reagent which indicates the endpoint or equivalent point of the titration. The

strength of concentration of a solution is expressed in the following ways.

NORMALITY: Number of gram equivalents of the substance dissolved per liter of the solution is called Normality. It is denoted by N

$$\text{Normality} = \frac{W_{\text{solute}}}{E_{\text{solute}}} \times \frac{1}{V_{\text{solvent}}} \quad (\text{in lit})$$

Where E is Gram equivalent weight

MOLARITY: Number of gram moles of a solute dissolved per liter of solution is called Molarity. It is denoted by M

$$\text{Molarity} = \frac{W_{\text{solute}}}{M_{\text{solute}}} \times \frac{1}{V_{\text{solvent}}} \quad (\text{in lit})$$

Where M is Gram molecular weight

MOLALITY: It is the number of moles of the substance dissolved in 1 kg of the solvent. It is denoted by (m).

$$\text{Molality} = \frac{W_{\text{solute}}}{M_{\text{solute}}} \times \frac{1}{W_{\text{solvent}}} \quad (\text{in kg})$$

Standard Solution:

A solution whose concentration is known is called a standard solution. Again there are two types of standard solutions depending on the nature of the substance.

1. Primary Standard Substances:

Any substance stable, pure, readily soluble in water, with high equivalent weight and the composition of its solutions should not change on standing or during storage is called a primary standard substance.

For example, crystalline oxalic acid, potassium dichromate and anhydrous sodium carbonate etc. are the primary standard substances.

2. Secondary Standard Substances:

The substance which does not fulfill the above mentioned requirements of primary standard substances, their solutions are not directly prepared by weighing and the exact strength of the solution is found by titrating it against some primary standard, called a secondary standard substance, and the process is called standardization.

Common secondary standard substances are NaOH, KOH which are hygroscopic, KMnO_4 which undergoes auto-decomposition in permanganate solution on standing and inorganic acids like HCl, H_2SO_4 whose concentrations are known approximately.

Types of Titrations

Depending on the nature of chemical reaction involved, the volumetric titrations are classified into the following types.

1. **Acid-Base Titrations**
2. **Oxidation-Reduction Titrations (Redox Titrations)**
3. **Complexometric Titrations**
4. **Precipitation titrations.**

1. **Acid-Base Titrations:** Acid-base titrations, in which an acidic or basic titrant reacts with an analyte that is a base or an acid.
2. **Oxidation-Reduction Titrations (Redox Titrations):** The titrations in which a reducing agent is titrated against an oxidizing agent and vice versa are called redox titrations.
3. **Complexometric Titrations:** The titrations involving the formation of a stable soluble complex between the metal and the complexing reagent. The metal ion is called the central atom and the complexing reagent is called the ligand.
4. **Precipitation titrations:** Precipitation titrations, in which the analyte and titrant react to form a precipitate.


DETERMINATION OF TOTAL HARDNESS OF WATER BY COMPLEXOMETRIC METHOD USING EDTA

Aim: To estimate the hardness of the given water sample by complexometric method using EDTA.


Instruments and Apparatus: Burette, pipette, conical flask, standard flask and burette stand.

Chemicals required: Magnesium Sulphate Heptahydrate ($\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$), Ethylene Diamine Tetra Acetic Acid (EDTA), Eriochrome Black T (EBT) indicator, Ammonia buffer solution, hard water and distilled water.


Principle: Hard water which contains Calcium and magnesium ions forms a wine red coloured complex with the Eriochrome Black-T indicator. EDTA forms a colourless stable complex with free metal ions like Ca^{2+} and Mg^{2+} .




When EDTA is added from the burette, it extracts the metal ions from the metal ion-indicator complex thereby releasing the free indicator (The stability of metal ion-indicator complex is less than that of the metal ion-EDTA complex, and hence EDTA extracts metal ion from the ion-indicator complex).



The reaction takes place at $\text{pH} = 10$ and the buffer is made by Ammonium chloride and Ammonia solution.



structure of EDTA



Structure of EBT

Procedure

Preparation of Standard $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ solution (0.05N):

Weight accurately 0.6162 g of Magnesium Sulphate and transfer it in to a clean 100 ml volumetric flask with the help of funnel, dissolved this crystal in minimum amount of distilled water and then made up to the mark with distilled water. It is shaken well to get uniform concentration and calculate the normality of the prepared Standard MgSO_4 solution.

Normality of Standard $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ solution = 0.05N

weight of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ = ?

246.48 = Molecular weight of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

123.24 = Equivalent weight of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

$$\text{Normality} = \frac{\text{Weight}}{\text{Equivalent weight}} \times \frac{1000}{\text{Volume (ml)}}$$

$$\text{Weight} = \text{Normality} \times \text{Gr. Eq. wt} \times \frac{\text{Volume (ml)}}{1000}$$

Titration I: Standardisation of EDTA Solution

Pipette out 20 ml magnesium sulphate solution in to a clean conical flask. To this add 2ml of buffer solution and 3 to 4 drops of EBT indicator. Then there develops a wine red colour. This solution is titrated against EDTA solution taken in the burette. The end point in this reaction is conversion of wine red colour into blue. Repeat the titrations till concurrent values are obtained. The normality of EDTA is calculated.

S.No.	Volume of Standard MgSO_4 (V_1) (ml)	Burette reading (ml)		Volume of EDTA Solution (V_2) (ml)
		Initial	Final	

Volume of standard MgSO_4 solution (V_1) = ml

Strength of standard MgSO_4 solution (N_1) = N

Volume of EDTA Solution (V_2) = ml

Strength of EDTA Solution (N_2) = N

$$N_1 V_1 =$$

$$N_2 V_2$$

$$N_1 V_1 / V_2$$

Pipette out 20 ml of hard water sample into a clean conical flask. To this add 2 ml of buffer solution followed by 3 to 4 drops of EBT indicator and titrate the solution till a clear blue colour persists. This marks the endpoint of the titration. Repeat the titrations for constant values and calculate the amount of total hardness present in the given sample.

S.No.	Volume of hard water (V ₃) (ml)	Burette reading (ml)		Volume of EDTA solution (V ₂) (ml)
		Initial	Final	

Strength of hard water solution (N₃) = N

Volume of EDTA Solution (V_2) = ml

Strength of EDTA Solution (N₂) = N

$$\mathbf{N}_2 \mathbf{V}_2 = \mathbf{N}_3 \mathbf{V}_3$$

$$N_3 = N_2 V_2 / V_3$$

Amount of total hardness present in the given sample = $N_3 \times 100 \times 1000 = \text{ppm}$

Transfer 100 ml of the given sample water into a beaker and boil it gently for 20 minutes. Cool and filter it directly into a 100 ml standard flask. Make up the solution with distilled water and mix well. Pipette out 20 ml of this solution into a clean conical flask. To this add 2 ml of buffer solution followed by 3 to 4 drops of EBT indicator and titrate the solution till a clear blue colour persists. This marks the endpoint of the titration. Repeat the titrations for constant values and calculate the amount of permanent hardness present in the given sample.

S.No.	Volume of sample water (V_4) (ml)	Burette reading (ml)		Volume of EDTA Solution (V_2) (ml)
		Initial	Final	

Volume of water solution (V_4) = ml

Strength of water solution (N_4) = N

Volume of EDTA Solution (V_2) = ml

Strength of EDTA Solution (N_2) = N

$$N_2 V_2 =$$

$$N$$

$$V_4 N_4 = N_2 V_2 / V_4$$

Calculation:

Amount of permanent hardness present in the given sample = $N_3 \times 100 \times 1000 = \text{ppm}$

Temporary hardness = Total hardness – Permanent hardness = ppm

Result: Total hardness of water sample = ppm

Permanent hardness of water sample = ppm

Temporary hardness of water sample = ppm

ESTIMATION OF ACID BY CONDUCTOMETRIC TITRATIONS


Aim: To determine the strength of given HCl solution by titrating against standard NaOH solution conductometrically.

Instruments and Apparatus: Burette, pipette, glass rod, beaker, standard flask, conductivity meter and conductivity cell with known cell constant.

Chemicals required: Hydrochloric acid (HCl), Oxalic acid and sodium hydroxide.

Principle:

According to Kohlrausch's law, electrical conductivity of a solution depends on the number of ions present in it. In the titration of strong acid HCl with strong base NaOH solution, before addition of NaOH solution there will be high conductance due to the presence of a large number of ions in it. Gradual addition of NaOH solution decreases the conductance due to the combination of H^+ ions with OH^- ions of the base to form undissociated water molecules. The conductance of the solution decreases till the neutralization point is reached and increases quickly due to free Na^+ and OH^- ions added in excess of NaOH solution. A plot of conductivity Vs volume of NaOH added will consist of two straight lines intersecting at the equivalent point. The shape of the plot is V.



Procedure:

Preparation of standard oxalic acid solution (0.1N):

Weigh 0.63 gram of oxalic acid crystals into a clean 100 ml standard flask, dissolve in a small amount of water and make up the solution up to the mark with distilled water. Shake the flask well for uniform concentration.

Normality of Oxalic Acid = 0.1N, W = weight of Oxalic Acid = ?

126.07 = Molecular weight of Oxalic Acid 63 = Equivalent weight of Oxalic Acid

$$\text{Normality} = \frac{\square\square\square\square\square - \square}{\square\square.\square\square.\square\square\square\square}$$

$$\text{Weight} = \text{Normality} \times \text{Gr. Eq. wt} \times \square - \square$$

Titration I: Standardisation of NaOH solution

Fill the burette with given NaOH solution. Pipette out 20 ml of oxalic acid into a conical flask and add 2 to 3 drops of phenolphthalein indicator and titrate the colourless solution against NaOH till pale pink colour is obtained as end point. Repeat the titration to get concurrent values.

S.No.	Volume of Oxalic acid solution (ml) (V_1)	Burette reading (ml)		Volume of NaOH (ml) (V_2)
		Initial	Final	

Volume of standard Oxalic acid solution (V_1) = _____ ml

Strength of standard Oxalic acid solution (N_1) =

N

Volume of Sodium hydroxide (V_2)

=

ml

Strength of Sodium hydroxide (N_2)

=

N

$$N_1 V_1$$

$$= N_2 V_2$$


$$N_2 = N_1 V_1 / V_2$$

Titration II: Determination of strength of HCl by conductivity meter

Transfer the given unknown HCl solution into a clean 100 ml standard flask and make up the solution upto mark with distilled water and mix well to obtain uniform concentration. Pipette out 20 ml of the given acid into a clean 100 ml beaker, dip the conductivity cell in it and measure the conductance initially. Fill the burette with standard NaOH solution and run down into the beaker (at a time 1 ml) with gentle stirring of the contents of the beaker and note the conductance after each addition. The measured conductance are recorded and tabulated. Take about 20 readings, by plotting the graph (conductivity Vs volume of NaOH) we can get the amount of NaOH required for the neutralization of HCl taken in the beaker (i.e. end point) and calculate the concentration of HCl solution from the known concentration of NaOH solution.

S.NO	Volume of NaOH added (ml)	Conductance (Ohm^{-1} or Siemen's)
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		

Model Graph:



Calculation:

Normality of NaOH (N_2) = N

Volume of NaOH (from graph) (V_2) = ml

Normality of HCl (N_3) = N

Volume of HCl (V_3) = ml

$$N_2 V_2 =$$

$$N$$

$$V_3 N_3 = N_2 V_2 / V_3$$

Strength of given HCl solution = Normality of HCl (N_3) \times Eq. Wt. of HCl (36.54)

Result: The strength of HCl = g / lit

Advantage of Conductometric Titrations:

1. Coloured solutions which cannot be titrated by ordinary volumetric titration can be carried out without the help of indicator.
2. This method can also be employed in case of very dilute solution and also for weak acid and weak base.
3. No special care is necessary near the endpoint as it is determined graphically.

ESTIMATION OF Fe+2 BY POTENTIOMETRY USING KMnO₄

Aim: To determine the amount of ferrous ion present in given whole solution by titrating against standard KMnO₄ solution potentiometrically.

Instruments and Apparatus: Burette, pipette, glass rod, beaker, standard flask, potentiometer, Platinum electrode and saturated calomel electrode (SCE).

Chemicals required: Ferrous Ammonium Sulphate (FAS), Potassium permanganate, Dil. Sulphuric acid.

Principle:

Redox titrations can be carried out potentiometrically using platinum-calomel electrode combination. For the reaction



The potential is given by Nernst equation,

$$E = E^\circ + \frac{RT}{nF} \log \frac{[\text{Oxidised form}]}{[\text{Reduction form}]}$$

Where, E° is the standard potential of the system. The potential of the system is thus controlled by the ratio of the concentration of the oxidized to that of the reduced species present. As the reaction proceeds, the ratio and hence the potential changes more rapidly in the vicinity of the endpoint of the titration. This may be followed potentiometrically and a plot of change in potential against volume (Titration curve) is characterized by sudden change in the potential at the equivalence point.

The reaction that takes place in the determination of Fe(II) is



The experiment cell to be used is:



Procedure:

Preparation of a Standard FAS solution (0.01N):

Weight accurately 0.392g of FAS and transfer it in to a clean 100 ml volumetric flask with the help of funnel, dissolve this crystal in minimum amount of distilled water and 2 ml of dil. Sulphuric acid and make up the solution up to the mark, shake well for uniform concentration and calculate the normality of the prepared FAS (Mohr's salt) solution.

Normality of FAS = 0.01, W = weight of FAS = ?

392.13 = Molecular weight of FAS

392.13 = Equivalent weight of FAS

$$\text{Normality} = \frac{\text{Weight of FAS}}{\text{Equivalent weight of FAS} \times \text{Volume of solution (ml)}}$$

$$\text{Weight} = \text{Normality} \times \text{Gr. Eq. wt} \times \text{Volume}$$

Titration I: Standardisation of KMnO₄

Pipette out 20 ml of standard FAS solution into a clean conical flask. Add one test tube full (about 20 ml) of dilute sulphuric acid solution. Titrate this mixture against potassium permanganate taken in the burette. The end point is the appearance of light pink violet colour. Repeat the titration for concordant titre values. From the strength of FAS calculate the strength of KMnO₄ solution.

S.No.	Volume of FAS (V ₁) (ml)	Burette reading (ml)		Volume of KMnO ₄ (V ₂) (ml)
		Initial	Final	

Volume of standard FAS solution (V₁) = ml

Strength of standard FAS solution (N₁) = N

Volume of KMnO₄ (V₂) = ml

Strength of KMnO₄ (N₂) = N

$$N_1 V_1 =$$

$$N$$

$$2 V_2 N_2 = N_1 V_1 / V_2$$


Calibration of Potentiometer:

Switch on the potentiometer and connect the standard cell terminals to either channel A (move channel switch to position A) or channel B (move the channel switch to position B). The meter should read 1.018 V. In case it is not 1.018 V, adjust the Std. knob to obtain reference value.

Titration II: Estimation of Ferrous ion


Transfer the given unknown ferrous ion solution into a clean 100 ml standard flask and make up the solution up to mark with distilled water and mix well to obtain uniform concentration. Pipette out 20 ml of above solution into a clean 100 ml beaker and add one test tube of Dil. Sulphuric acid. Place the electrode assembly (platinum electrode as indicator electrode and a SCE as reference electrode) in the beaker and connect the same to the potentiometer. Add permanganate from burette in 1 ml portions to the ferrous solution, stir it and note down the EMF values. Continue the titration till a sudden inflexion in EMF occurs, then take about 6 to 8 readings after in 1 ml intervals. From the titration approximate volume of permanganate required is found out. Draw a graph of E cell Vs volume of permanganate added, the intersecting point gives an approximate equivalence point.

S.No	Volume of KMnO_4 (ml)	EMF (volt)



The titration is repeated with addition of permanganate in 0.2 ml lots in the vicinity of end point (in 2 ml range). Plot a graph $\Delta E / \Delta V$ Vs volume of KMnO_4 added. Calculate the normality of ferrous solution and determine the amount of iron in the given solution.

S.No	Volume of KMnO ₄ (ml)	EMF(volt)	ΔE (volt)	ΔV (ml)	$\Delta E/\Delta V$ (volt/ml)	Average Volume (ml)



Calculation:

Normality of KMnO_4 (N_2) = N

Volume of KMnO_4 (from graph) (V_2) = ml

Normality of FAS (N_3) = N

Volume of FAS (V_3) = ml

$$N_2 V_2 =$$

$$N$$

$$V_3 N_3 = N_2 V_2 / V_3$$

Weight of FAS per litre = Normality of FAS (N_3) \times Eq. Wt. of FAS (55.85)

Weight of FAS per 100 ml = Weight of FAS per litre / 10

Result: The weight of Fe(II) present in the given solution by potentiometry is found to be _____ gr/lit

DETERMINATION OF STRENGTH OF GIVEN HCl USING pH METER

AIM


To determine the strength of given HCl using pH meter, a standard solution of NaOH of N is provided.

PRINCIPLE

Since the pH of the solution is related to the H⁺ ion concentration by the following formula,

$$\text{pH} = -\log [\text{H}^+]$$

measurement of pH of the solution gives the concentration of H⁺ ions in the solution. When NaOH is added slowly from the burette to the solution of HCl, the fast moving H⁺ ions are progressively replaced by slow moving Na⁺ ions. As a result pH of the solution increases.




The increase in pH takes place until all the H⁺ ions are completely neutralised (upto the end point). After the end point, further addition of NaOH increases the pH sharply as there is an excess of fast moving OH⁻ ions.

MATERIALS REQUIRED


(i) pH meter, (ii) Glass electrode, (iii) 100 Beaker, (iv) Standard NaOH, (v) Given HCl, (vi) Burette, pipette, Glass rod etc., (vii) Distilled water.

S.No.	Volume of NaOH (ml)	pH	ΔpH	ΔV (ml)	$\Delta pH / \Delta V$
1.					
2.					
3.					
4.					
5.					
6.					
7.					
8.					
9.					
10.					
11.					
12.					
13.					
14.					
15.					
16.					
17.					
18.					
19.					
20.					

Graph I (Titration I)



Graph II (Titration II)



Calculation

Step I: Calculation of Strength of HCl

Volume of HCl, $V_1 = 20$ ml

Strength of the HCl, N_1 = ?

Volume of the NaOH, V_2 = (titre value)

Strength of the NaOH, N_2 = N

According to the law of volumetric analysis, $V_1N_1 = V_2N_2$

$$N_1 = \frac{V_2 N_2}{V_1}$$
$$= \frac{\text{.....} \times \text{.....} N}{20}$$

Strength of HCl present in 1 litre of the given solution = N

Step II: Calculation of amount of HCl

$$\left. \begin{array}{l} \text{The amount of HCl present in} \\ \text{1000 ml of the given solution} \end{array} \right\} = \text{.....} N \times \text{Eq. wt. of HCl (36.45)}$$
$$= \text{.....} 36.45$$
$$= \text{.....} \text{ gms}$$

Step II: Calculation of amount of HCl

The amount of HCl present in 1000 ml of the given solution = N \times Eq. wt. of HCl (36.45)

= 36.45

= gms

PROCEDURE

TITRATION - I


The given hydrochloric acid solution is transferred into 100 ml standard flask and made up to the zero mark using distilled water. 20 ml of this made up solution is pipetted out into a clean 100 ml beaker and diluted by adding 20 ml of distilled water. The glass electrode is dipped in it and connected with a pH meter.

The burette is washed well with water and rinsed with a small amount of given NaOH solution. It is then filled with the same upto zero mark.


Titration is carried out by adding std. NaOH solution in portions of 1 ml from the burette to the HCl solution taken in the beaker and pH of the solution is noted for each addition. This process is continued until at least 5 readings are taken after the end point, and the range at which the end point lies is found out by plotting volume of NaOH added against pH (graph I).

TITRATION - II

Another titration is carried out by adding std. NaOH solution in portions of 0.1 ml near the end point and pH of the solution is noted after each addition. The addition of NaOH is continued even after the end point for further 1 ml. The accurate end point is found out by plotting $\Delta\text{pH} / \Delta V$ against volume of NaOH added (graph II). From the end point, the strength of HCl solution and hence the amount of HCl is calculated.



Graph :



RESULT

1. Strength of the given HCl solution = N
2. Amount of HCl present in 1 litre of the solution = gms.

PREPARATION OF BAKELITE

AIM: - To prepare phenol formaldehyde resin. (Bakelite)


CHEMICALS USED: -

Glacial acetic acid or 2N NaOH,
40% formaldehyde solution,
Phenol,
conc. H₂SO₄ /HCl

APPARATUS REQUIRED: -Glass rod, beakers, funnel, measuring cylinder, dropper and filter paper.

PRINCIPLE: - Phenol formaldehyde resins (PFs) are condensation polymers and are obtained by condensing phenol with formaldehyde in the presence of an acidic or alkaline catalyst. They were first prepared by **Backeland**, an American Chemist who gave them the name as **Bakelite**. These are thermosetting polymers.

PREPARATION:- PFs are prepared by reaction of phenol with formaldehyde in the presence of acidic or basic catalyst. The process may be carried out as follows
A mixture of phenol and formaldehyde are allowed to react in the presence of a catalyst. The process involves formation of methylene bridges in *ortho*, *para* or *both ortho and para* positions. This results first in the formation of linear polymer (Called **NOVALAC**) and then in to cross-linked polymer called phenol-formaldehyde resin or **Bakelite**



PROCEDURE: -

1. Place 5ml of glacial acetic acid and 2.5ml of 40% formaldehyde solution in a 500ml beaker and add 2 grams of phenol.
2. Add few ml of conc. Sulphuric acid into the mixture carefully. Within 5 min. a large mass of plastic is formed.
3. The residue obtained is washed several times with distilled water, and filtered product is dried and yield is calculated.

USES:-

They are used for making moulded articles such as radio and TV parts, combs, fountain pen barrels, phonograph records etc.

PRECAUTIONS:

1. The reaction is sometimes vigorous and it is better to be a few feet away from the beaker while adding the H_2SO_4 and until the reaction is complete.
2. The experiment should be preferably carried out in fume cupboard.

RESULT: - The weight of the phenol formaldehyde resin is ____ g

Preparation Nylon – 6

AIM: -

To synthesize Nylon-6 polymer through ring-opening polymerization of caprolactam.


CHEMICALS USED: -

- ϵ -Caprolactam (monomer)
- · Acetic acid (to terminate the reaction)
- · Heat source (oil bath or hot plate)
- · Glass rod
- · Beaker or test tube
- · Ice bath
- · Distilled water

PRINCIPLE: -

Nylon-6 is a **synthetic polyamide** made by the **ring-opening polymerization** of ϵ -caprolactam. Unlike Nylon-6,6, which is made from two monomers (hexamethylenediamine and adipic acid), Nylon-6 is derived from a single monomer. This reaction is typically catalyzed by heat.

Reaction:



Procedure:

Take 5–10 grams of ϵ -caprolactam in a clean, dry test tube or small beaker.

Place the test tube in an oil bath and heat it to around 220–250°C. Stir gently using a glass rod.

Maintain this temperature for about 15–20 minutes until a viscous melt is formed.

Cool the melt slowly at room temperature and then place it in an ice bath.

The solid polymer that forms is Nylon-6.

Wash the solid polymer with cold water and a few drops of acetic acid to remove unreacted monomer and neutralize any residual base.

Dry the product and note its appearance and texture.

RESULT: - The weight of the Nylon-6 is ____ g.

Precautions:

Handle hot apparatus with care.

Perform the experiment in a well-ventilated area or under a fume hood.

Use gloves and goggles as caprolactam may irritate skin.

DETERMINATION OF ACID VALUE OF COCONUT OIL

Aim: To determine the acid value present in given Coconut oil.

Chemicals: Oxalic acid, KOH Solution, Alcohol, Coconut oil, Phenolphthalein indicator and distilled water.

Apparatus: Burette, pipette, conical flask, Wash bottle, Dropper.

Principle: Acid value indicates the proportion of free fatty acid present in an oil or fat and may be defined as the number of milligrams of caustic potash required to neutralize the acid in 1 gm of the sample. The normal acid value for most samples lies within 0.5. If any titratable acid other than a fatty acid is present in the sample, it will be an error. A high acid value indicates a stale oil or fat stored under improper conditions.

Procedure

Standardization of KOH: Pipette out 20 ml of 0.1 N oxalic acid solution in a 250 ml conical flask. Add 1 or 2 drops of phenolphthalein indicator to this solution. Titrate this solution against KOH taken in a burette. The appearance of pink color indicates the end point. From the volume of the KOH solution in burette, find the normality of KOH.

Determination of Acid value in coconut oil: Weigh 5 gm of coconut oil and transfer it into 250 ml conical flask. Add 50 ml of neutralized alcohol solution to the oil solution. Heat this mixture for 10 minutes by using the heater. Take the solution after 10 minutes and add 1 or 2 drops of phenolphthalein indicator. Titrate this against the KOH solution from the burette. The appearance of pink color indicates the end point.

Calculation

Titration I: Standardisation of Potassium hydroxide solution:

S.No	Volume of Oxalic acid (ml)	Burette Reading		Concordant volume of KOH (ml)
		Initial	Final	

Burette: KOH

Pipette solution: Oxalic acid Indicator:

Phenolphthalein

Endpoint: Appearance of pink colour

Volume of oxalic acid (V_1) = Normality of

oxalic acid (N_1) =

Volume of KOH consumed (V_2) =

Normality of KOH consumed (

$N_2 = V_1 N_1 / V_2$ Normality of KOH (N_2) = ____.

Titration II:

Estimation of acid value

S.No	Volume of Oxalic acid (ml)	Burette Reading		Concordant Volume of KOH (ml)
		Initial	Final	

Burette solution: KOH

Pipette solution: Oil + 50 ml of neutralized alcohol

Indicator: Phenolphthalein

Endpoint: Appearance of pink colour

Acid value = $(\text{Volume of KOH} \times \text{Normality of KOH} \times \text{Eq. wt} \times 1000) / \text{Weight of Oil sample}$

Acid Value = _____.

RESULT

The acid value of the given Coconut oil sample was found to be _____.

DETERMINATION OF VISCOSITY OF CASTOROIL AND GROUNDNUT OIL BY USING OSTWALD'S VISCOMETER.


Aim: To determine the viscosity of a given liquid sample by Oswald's viscometer.

Instruments and Apparatus: Oswald's viscometer and stopwatch.

Chemicals required: Liquid sample and water sample.

Principle: Viscosity is the most important property of any lubricating oil. It describes the suitability of liquid for lubrication purposes. Viscosity of any liquid is that characteristic of fluid by virtue of which it hinders its own flow. Absolute viscosity is a measure of the internal resistance. Absolute viscosity is the tangential force per unit area required to move one horizontal plane with respect to the other at unit velocity when maintained at a unit distance apart by the fluid. Kinematic viscosity is the ratio of absolute or dynamic viscosity to density - a quantity in which no force is involved. Kinematic viscosity can be obtained by dividing the absolute viscosity of a fluid with its mass density.

Viscosity index (vi) is an arbitrary measure for the change of viscosity with temperature. It is used to characterize lubricating oil in the automotive industry.



Procedure:

Determination of density of given sample

Weight of specific gravity bottle + given sample $W_1 =$ g

Weight of specific gravity bottle $W_2 =$ g

Weight of given sample $(W_1 - W_2) =$ g

Volume of given sample $(V) =$ ml

Density of given sample $(d_2) = W_1 - W_2 / V (\text{gm/cc}) =$

Take a clean and dryOswald viscometer and set it vertically on a stand. Introduce 10 ml of water throughE(refer figure) intothe largebulb(B). Sucktheliquidup intothebulb‘A’througharubber tubing attached to the end ‘D’to a levelabove the mark ‘X’. Allow the liquid to flow freelythrough thecapillaryand notethetime t_1 (usethewrist watchor stop watch) fortheliquid flow fromXto Y. repeat stepstwoormoretimesandgetanaveragevalueoft₁.Repeattheprocесstoobtainthe averagetime for the liquids to flow fromX to Y. Tabulate the data as follows.

S.No.	Sample	Density	Timeflowtime average		Timeaverage
			Exp1	Exp2	
1	Water	1			
2					
3					

Calculations:

$$\eta_1 / \eta_2 = (d_1 t_1) / (d_2 t_2)$$

Where η_1 = viscosity of water,

η_2 = viscosity of given liquid d_1 =

density of water,

d_2 = density of given liquid

t_1 = time flow of water from X to Y

t_2 = time flow of given liquid from X to Y

Result: The viscosity of the given liquid = (unit of viscosity are poise, Pascal - second or stokes)

Determination of Rate of Corrosion of Mild Steel in the Presence and Absence of Inhibitor

Aim:

- Mild steel specimens (pre-weighed and cleaned)
- 1N HCl solution (as the corrosive medium)
- Corrosion inhibitor (e.g., sodium benzoate, hexamine, or potassium dichromate)
- Beakers
- Digital balance (accuracy: 0.001 g)
- Stop watch or timer
- Glass rod
- Tongs and gloves

Principle: Corrosion is the electrochemical deterioration of metals due to reactions with their environment. The corrosion rate can be calculated using the weight loss method:

Corrosion Rate (mg/cm²/day) = $\frac{\text{Area} \times \text{Time} (\text{Initial Weight} - \text{Final Weight})}{\text{Final Weight}}$

Corrosion inhibitors reduce the rate of corrosion by forming a **protective film** on the metal surface or altering the electrochemical reaction.

Procedure:

Step 1: Surface Preparation

Clean the mild steel specimens with emery paper.

Wash with distilled water, then with acetone or alcohol, and dry.

Measure and record the initial weight (W_1) using a digital balance.

Measure the surface area of each specimen.

Step 2: Without Inhibitor

Place one specimen in a beaker containing 100 mL of 1N HCl.

Leave it undisturbed for a fixed duration (e.g., 60 minutes).

Step 3: With Inhibitor

Prepare a second beaker with **100 mL of 1N HCl + known concentration of inhibitor**.

Immerse the second specimen for the **same duration**.

Step 4: Post-Exposure

After the exposure time, remove both specimens, clean off corrosion products using a soft brush, rinse with water, dry, and **reweigh (W_2)**.

Calculate the **weight loss ($W_1 - W_2$)** and use it to determine the corrosion rate.

Calculations:

weight Loss (mg) = $W_1 - W_2$

Corrosion Rate = $\frac{\text{Weight Loss (mg)}}{\text{Surface Area} \times \text{Time (days)}}$

Compare the corrosion rates **with and without the inhibitor**.

Result:

The corrosion rate without inhibitor was found to be: ____ mg/cm²/day

The corrosion rate with inhibitor was found to be: ____ mg/cm²/day

The inhibitor **effectively reduced** the corrosion rate of mild steel.

Virtual lab experiments:

Construction of Fuel cell and its working.

<https://youtu.be/kIPMhijlchs>

Smart materials for Biomedical applications

https://youtu.be/QYkjFTgF8_k

Batteries for electrical vehicles

<https://youtu.be/ZIJMwyLfO88>

Functioning of solar cell and its applications

<https://youtu.be/IODbzD8LiKM>